


A1



A2



A3



A4



A5



A6



A7



A8



A8-2



A9



A10



A11



A12



A13

### 产品细节



L-A



L-B



L-C

品名	型号	适用模块
透明盖	L-A	A1-A9
透明盖	L-B	A10-A11
透明盖	L-C	A12-A14

\* 单管、排管配 L-A；带裙边 PCR 板配 L-B

## GH/GC-100 系列可更换模块

### 可更换模块

型号	孔径及试管数	最高温度	型号	孔径及试管数	最高温度	型号	孔径及试管数	最高温度
A1	96×0.2ml 标准板	100°C	A6	24×φ12mm 试管	100°C	A10	96×0.2ml 酶标板 (平底)	100°C
A2	54×0.5ml 离心管	100°C	A7	32×0.2ml+25×1.5ml 离心管	100°C	A11	24×5ml	100°C
A3	35×1.5ml 离心管	100°C	A8	32×0.2ml+10×0.5ml+15×1.5ml 离心管	100°C	A12	12×15ml	100°C
A4	35×2.0ml 离心管	100°C	A8-2	32×0.2ml+10×0.5ml+15×2.0ml 离心管	100°C	A13	6×50ml	100°C
A5	20×0.5ml+15×1.5ml 离心管	100°C	A9	103×67×30 (方槽模块)	100°C	A14	12×10ml	100°C