

GA150 系列可更换模块

/// 温控系列



///D1



///D2



///D3



///D4



///D5



///D6



///D7



///D8




///D9



///D10



///D11



///D12



///D13




///D14



///D15



///D16



///D17



///D18

/// 可更换模块

型号	试管直径	试管数	模块尺寸 (W×D×H)	型号	试管直径	试管数	模块尺寸 (W×D×H)
D1	6mm	48	96.5×76.5×50 (mm)	D10	26mm	8	96.5×76.5×50 (mm)
D2	7mm	48	96.5×76.5×50 (mm)	D11	28mm	6	96.5×76.5×50 (mm)
D3	10mm	24	96.5×76.5×50 (mm)	D12	38mm	3	96.5×76.5×50 (mm)
D4	12mm	24	96.5×76.5×50 (mm)	D13	0.2ml 离心管	48	96.5×76.5×50 (mm)
D5	13mm	24	96.5×76.5×50 (mm)	D14	0.5ml 离心管	48	96.5×76.5×50 (mm)
D6	15mm	12	96.5×76.5×50 (mm)	D15	1.5ml 离心管	24	96.5×76.5×50 (mm)
D7	16mm	12	96.5×76.5×50 (mm)	D16	2.0ml 离心管	24	96.5×76.5×50 (mm)
D8	19mm	12	96.5×76.5×50 (mm)	D17	0.2ml	96	96.5×153×33.5 (mm)
D9	20mm	12	96.5×76.5×50 (mm)	D18	平板	酶标板、深孔板、细胞培养板	96.5×153×22.5 (mm)